



To verify that this document is original, please check this QR code

VOL-01

PERFORMANCE QUALIFICATION PROTOCOL

Annex-01.01

Results Data Sheets

FOR
AUTOCLAVE
(M-AUX-AU014)
OF
MIDDLE EAST FOR VACCINES (MEVAC)
AL-SALEHIA

CLIENT:

MEVAC

QUALIFICATION PROVIDER:



Package Content:

VOL-01 Performance Qualification Protocol

Annex-01.01) Results Data Sheets

Annex-01.02) Turn Over Package (Attachments)

EGYPT: 26 Al Israa St., Lebanon Sq., Mohandseen, Giza, Egypt

KSA: 3321, Prince Mamdouh Bin Abdulaziz St., Al Sulaymaniyah, Riyadh, KSA

UAE: Sharjah Media City, Sharjah, UAE

www.technoairgate.com

info@technoairgate.com

+201023310801

+966543728852

+971502013049

+20233034440



Table of Content

	Page
1. REPORT REVIEW AND APPROVAL	3
2. INTRODUCCION	4
2.1. Report Purpose	4
2.2. System Brief Description	4
2.3. Associated Material and Equipment Description	4
2.4. Study Design	5
3. RESULTS DATA SHEETS	6
DATA SHEET (01.01.01) MAIN CYCLES TEST RESULTS	6
01.01.01.01) MAXIMUM SOLID DISTRIBUTION CYCLE TEST RESULTS	7
01.01.01.01.01) Maximum solid distribution Cycle Test Results Data Sheets	8
01.01.01.01.02) Maximum solid distribution Cycle Conclusion and Recommendation	9
01.01.01.01.03) Maximum solid distribution Cycle Layout and Sensors' Distribution	10
01.01.01.01.04) Maximum solid distribution Cycle Live Photos	11
01.01.02.01) MAXIMUM SOLID PENETRATION CYCLE TEST RESULTS	12
01.01.01.02.01) Maximum solid penetration Cycle Test Results Data Sheets	13
01.01.01.02.02) Maximum solid penetration Cycle Conclusion and Recommendation	14
01.01.01.02.03) Maximum solid penetration Cycle Layout and Sensors' Distribution	15
01.01.01.02.04) Maximum solid penetration Cycle Live Photos	16
DATA SHEET (01.01.02) DEVIATION REPORT	17
01.01.02.01) Deviation Report	18
DATA SHEET (01.01.03) REFERENCE EQUIPMENT LIST AND CALIBRATION VERIFICATION	19
01.01.03.01) Reference Equipment List	20
01.01.03.02) Sensors Pre-Calibration Check	21
01.01.03.03) Sensors Post-Calibration Check	22
DATA SHEET (01.01.04) LIST OF ATTACHEMENTS	23
01.01.04.01) List of Attachments	24
Final Review And Approval	25



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 3 of 25

1. REPORT REVIEW AND APPROVAL

TAG Author(s)

Name	Title	Company	Date	Signature
Salma Salah	T.O.E.	TAG	30-06-2025	Salma Salah

TAG Inspector(s)

Name	Title	Company	Date	Signature
Ahmed Mohamed	P.E.	TAG	30-06-2025	Ahmed

TAG Reviewers(s)

Name	Title	Company	Date	Signature
Noha Essam	SR. T.O.E	TAG	30-06-2025	Noha Essam

TAG Approver(s)

Name	Title	Company	Date	Signature
Ahmed Tarek	T.M.	TAG	30-06-2025	Ahmed

MIDDLE EAST FOR VACCINES (MEVAC) Reviewer(s)

Name	Title	Company	Date	Signature

MIDDLE EAST FOR VACCINES (MEVAC) Approver(s)

Name	Title	Company	Date	Signature

Change History

No.	Revision No.	Amended Page No.	Reason for change	Issued Date
1.	00	NA	Original	30-06-2025

Note: Any amendment will be indicated with italic and bold font.



2. INTRODUCCION

2.1. Report Purpose

The purpose of this report is to document the results of the performance qualification study conducted to evaluate the heat distribution within the Autoclave located at MIDDLE EAST FOR VACCINES (MEVAC). The report includes the following points:

2.1.1 Detailed tables for the main cycle(s) conducted presenting the recorded:

- Average, maximum, and minimum recorded temperature for each location.
- F₀-value for each location.
- Average temperature across all locations.

2.1.2 Identifying the hottest and coolest points and illustrating their locations on the layout diagram.

2.1.3 Identifying any excursions or abnormal variations observed during the study and providing a summary of those deviations in the deviation report along with any corrective action taken.

2.2. System Brief Description

The scope of this report intends to cover the performance qualification results of Autoclave (M-AUX-AU014) located at MIDDLE EAST FOR VACCINES (MEVAC) with the following description:

System Brief Description	
Description	Autoclave
Dimensions/Volume	53 Liter
Manufacturing	TOMY
Model	SX-500
S.N.	53133024
Code	M-AUX-AU014
Location	Washing B.1.2.240

2.3. Associated Material and Equipment Description

Used Reference Equipment Description						
Description	No. of sensors	Range	Accuracy	Manufacturer	Model	Calibration Due Date
Data Loggers	10 sensors	(-40) °C to (140) °C	± 0.5 °C	Tecnosoft s.r.l.	SterilDisk Probe	Nov., 2025
Data Loggers	01 sensors	(0) bar to (4) bar	± 5% of F.S	Tecnosoft s.r.l.	PressureDisk 05	Nov., 2025



2.4. Study Design

Heat Distribution Study Design

Sterilization Temp. Setting	122 °C
Accepted Sterilization Temp.	(121 - 125) °C
Sterilization Time Setting	for 20 min
Accepted Sterilization Time	NLT 20 min
Accepted F ₀ -value	≥ 15 (as per customer requirement)

Cycle No.	Cycle Description	Load Description	Start Date	Stop Date	Sample Frequency
1.	Max. Solid distribution	2 Bottles (200 ml) & silicon tools	26/06/2025 12:26:00 PM	26/06/2025 02:08:15 PM	15 seconds
2.	Max. Solid penetration	2 Bottles (200 ml) & silicon tools with BI	26/06/2025 02:25:00 PM	26/06/2025 04:08:15 PM	15 seconds



3. RESULTS DATA SHEETS**VOL-01) PERFORMANCE QUALIFICATION PROTOCOL****ANNEX 01.01****RESULTS DATA SHEETS****DATA SHEET 01.01.01****MAIN CYCLES TEST RESULTS****Content:**

01.01.01.01	Maximum Solid Distribution Cycle Test Results
01.01.02.01	Maximum Solid Penetration Cycle Test Results



VOL-01) PERFORMANCE QUALIFICATION PROTOCOL

ANNEX 01.01

RESULTS DATA SHEETS

DATA SHEET 01.01.01.01

SOLID DISTRIBUTION CYCLE TEST RESULTS



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 8 of 25

Datasheet number 01.01.01.01.01: Maximum Solid distribution Cycle Test Results Data Sheets

Heat Distribution Study Data Sheet for Autoclave (M-AUX-AU014)

Cycle No.	1 of 2		
Cycle Name	Maximum Solid distribution Cycle		
Load Description	2 Bottles (200 ml) & silicon tools		
Sterilization Temp. Setting	122 °C		
Accepted Sterilization Temp.	(121 - 125) °C		
Sterilization Time Setting	for 20 min		
Accepted Sterilization Time	NLT 20 min		
Accepted F ₀ -value	≥ 15 (as per customer requirement)		
Cycle Start Time	26/06/2025 12:26:00 PM	Cycle End Time	26/06/2025 02:08:00 PM
Sterilization Start Time	26/06/2025 01:14:00 PM	Sterilization End Time	26/06/2025 01:34:30 PM

Sensor Position	Sensor S/N	Actual			F ₀ -value	Within Acceptance Criteria (YES / NO)?
		Average	Maximum	Minimum		
		Temp. (°C)	T (°C)	T (°C)		
1.	1011004332	122.74	122.90	121.73	30.33	YES
2.	1001009966	122.06	122.21	121.07	25.92	YES
3.	1011004203	122.68	122.81	121.63	29.90	YES
4.	1011004509	122.67	122.81	121.61	29.83	YES
5.	1011004211	122.41	122.53	121.42	28.10	YES
6.	1011004106	122.75*	122.91*	121.71	30.36	YES
7.	1011004102	122.60	122.71	121.56	29.31	YES
8.	1011004404	122.72	122.85	121.66	30.17	YES
9.	1011004338	121.90**	122.03	121.03**	24.97**	YES
10.	1011004503	122.63	122.77	121.61	29.52	YES
11.	Pressure	2.22	2.23	2.12	--	--
Average for Autoclave		122.52				

*: Maximum & Hottest.

** : Minimum & Coolest.

NOTES:

Reviewed By (TAG)

Reviewed By (MEVAC)

Name

Ahmed Tarek

Signature



Date

30-06-2025



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 9 of 25

Datasheet number 01.01.01.01.02: Maximum Solid distribution Cycle Final Report and Conclusion

- All the registered charts, graphs, and raw data (CD) are attached in the TOP.
- The following table presents a summary of the monitored data all over the sterilization period. All obtained values are not rounded.

Results Summary			
Parameter	Value	Location No.	Sensor S.N.
Min. Sterilization Temp.	121.03 °C	9	1011004338
Max. Sterilization Temp.	122.91 °C	6	1011004106
Hottest Point Temperature	122.75 °C	6	1011004106
Coolest Point Temperature	121.90 °C	9	1011004338
Minimum F0-value	24.97	9	1011004338
Slowest Point		1011004338	
Equilibrium Time		45 sec	
Holding Time (Exposure Time)		20 min & 30 sec	
Plateau Period		21 min & 15 sec	

Reviewed By (TAG)

Reviewed By MEVAC

Name

Ahmed Tarek


Signature



Date

30-06-2025



**Datasheet number 01.01.01.01.03: Maximum Solid distribution Cycle Layout and Sensors'
Distribution**

*: Hottest.
**: Coolest.



Datasheet number 01.01.01.01.04: Maximum Solid distribution Cycle Live Photos



VOL-01) PERFORMANCE QUALIFICATION PROTOCOL**ANNEX 01.01****RESULTS DATA SHEETS****DATA SHEET 01.01.01.01****SOLID PENETRATION CYCLE TEST RESULTS**

Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 13 of 25

Datasheet number 01.01.01.02.01: Maximum Solid Penetration Cycle Test Results Data Sheets

Heat Distribution Study Data Sheet for Autoclave (M-AUX-AU014)

Cycle No.	2 of 2
Cycle Name	Maximum Solid Penetration Cycle
Load Description	2 Bottles (200 ml) & silicon tools with BI
Sterilization Temp. Setting	122 °C
Accepted Sterilization Temp.	(121 - 125) °C
Sterilization Time Setting	for 20 min
Accepted Sterilization Time	NLT 20 min
Accepted F ₀ -value	≥ 15 (as per customer requirement)

Cycle Start Time	26/06/2025 02:25:00 PM	Cycle End Time	26/06/2025 04:08:15 PM
Sterilization Start Time	26/06/2025 03:03:45 PM	Sterilization End Time	26/06/2025 03:24:15 PM

Sensor Position	Sensor S/N	Actual			F ₀ -value	Within Acceptance Criteria (YES / NO)?
		Average	Maximum	Minimum		
		Temp. (°C)	T (°C)	T (°C)		
1.	1011004332	122.58	122.67	121.76	29.18	YES
2.	1001009966	122.04	122.12	121.30	25.77	YES
3.	1011004203	122.53	122.62	121.71	28.87	YES
4.	1011004509	122.53	122.62	121.70	28.86	YES
5.	1011004211	122.33	122.40	121.64	27.56	YES
6.	1011004106	122.70*	122.80*	121.93	30.02	YES
7.	1011004102	122.55	122.63	121.77	28.98	YES
8.	1011004404	122.59	122.69	121.81	29.25	YES
9.	1011004338	121.80**	121.86	121.17**	24.36**	YES
10.	1011004503	122.57	122.69	121.80	29.12	YES
11.	Pressure	2.21	2.22	2.15	--	--
Average for Autoclave		122.42				

*: Maximum & Hottest.

** : Minimum & Coolest.

NOTES:



Reviewed By (TAG)

Reviewed By (MEVAC)

Name

Ahmed Tarek

Signature

Date

30-06-2025

Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 14 of 25

Datasheet number 01.01.01.02.02: Maximum Solid Penetration Cycle Final Report and Conclusion

- All the registered charts, graphs, and raw data (CD) are attached in the TOP.
- The following table presents a summary of the monitored data all over the sterilization period. All obtained values are not rounded.

Results Summary			
Parameter	Value	Location No.	Sensor S.N.
Min. Sterilization Temp.	121.17 °C	9	1011004338
Max. Sterilization Temp.	122.80 °C	6	1011004106
Hottest Point Temperature	122.70 °C	6	1011004106
Coolest Point Temperature	121.80 °C	9	1011004338
Minimum F0-value	24.36	9	1011004338
Slowest Point		1011004338	
Equilibrium Time		30 sec	
Holding Time (Exposure Time)		20 min & 30 sec	
Plateau Period		21 min	

Reviewed By (TAG)

Reviewed By (MEVAC)

Name

Ahmed Tarek

Signature




Date

30-06-2025



Datasheet number 01.01.01.02.03: Maximum Solid Penetration Cycle Layout and Sensors' Distribution



*: Hottest.
**: Coolest.



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date:02/02/2025

Page 16 of 25

Datasheet number 01.01.01.02.04: Maximum Solid Penetration Cycle Live Photos



VOL-01) PERFORMANCE QUALIFICATION PROTOCOL

ANNEX 01.01


RESULTs DATA SHEETs

DATA SHEET 01.01.02

DEVIATION REPORT



Datasheet Number 01.01.02.01: Deviation Report

DEVIATION REPORT		Copy 01 of 01		
Test Report (Nº. and Description):		Deviation Nº.:		
Description:				
				
			Signature TAG:	Date:
			Signature MIDDLE EAST FOR VACCINES (MEVAC):	Date:
			Corrective Action:	
Person in charge to close the deviation		Expected closure date:		
Signature TAG:		Date:		
Signature QA MIDDLE EAST FOR VACCINES (MEVAC):		Date:		
Corrective Action results				
Signature MIDDLE EAST FOR VACCINES (MEVAC):		Date:		
Deviation successfully solved?	YES _____	NO _____		
Signature:		Date:		
Signature:		Date:		



VOL-01) PERFORMANCE QUALIFICATION PROTOCOL

ANNEX 01.01

RESULTS DATA SHEETS

DATA SHEET 01.01.03

REFERENCE EQUIPMENT LIST AND CALIBRATION VERIFICATION



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
 Issue Date: 02/02/2025

Page 20 of 25

Datasheet Number 01.01.03.01: Reference Equipment List

Reference Equipment List					
Instrument Description				CALIBRATION	
No.	Instrument Description	Manufacturer/ Model	S.N.	Most recent calibration date	Calibration Due Date
1.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004332	May., 2025	Nov., 2025
2.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1001009966	May., 2025	Nov., 2025
3.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004203	May., 2025	Nov., 2025
4.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004509	May., 2025	Nov., 2025
5.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004211	May., 2025	Nov., 2025
6.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004106	May., 2025	Nov., 2025
7.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004102	May., 2025	Nov., 2025
8.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004404	May., 2025	Nov., 2025
9.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004338	May., 2025	Nov., 2025
10.	Data Logger	Tecnsoft s.r.l./ SterilDisk Probe	1011004503	May., 2025	Nov., 2025
11.	Data Logger	Tecnsoft s.r.l./ PressureDisk 05	3501009327	May., 2025	Nov., 2025

Reviewed By (TAG)

Name

Ahmed Tarek

Signature

Date

30-06-2025

Reviewed By (MEVAC)



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 21 of 25

Datasheet Number 01.01.03.02: Sensors Pre-Calibration Check

Reference Equipment List

Equipment Description	Manufacturer	Model	S.N.	Most recent calibration date	Calibration Due Date
Pt-100 with digital indicator	GREISINGER	GMH 3710	32403245/116	Mar., 2025	Mar., 2026

Pre-Calibration Check Results

No.	Sensors S.N.	Error @ 100 °C	Error @ 121 °C	Error @ 130 °C
		(°C)	(°C)	(°C)
1.	1011004332	-0.08	0.11	0.07
2.	1001009966	0.13	0.26	0.22
3.	1011004203	0.19	0.13	0.05
4.	1011004509	0.39	0.19	0.22
5.	1011004211	0.46	0.48	0.43
6.	1011004106	-0.01	0.14	0.11
7.	1011004102	-0.04	0.14	0.09
8.	1011004404	0.45	0.24	0.29
9.	1011004338	0.03	0.16	0.27
10.	1011004503	0.15	0.26	0.20

Reviewed By (TAG)

Reviewed By (MEVAC)

Name

Ahmed Tarek

Signature

Date

30-06-2025



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 22 of 25

Datasheet Number 01.01.03.03: Sensors Post-Calibration Check

Reference Equipment List

Equipment Description	Manufacturer	Model	S.N.	Most recent calibration date	Calibration Due Date
Pt-100 with digital indicator	GREISINGER	GMH 3710	32403245/116	Mar., 2025	Mar., 2026

Post-Calibration Check Results

No.	Sensors S.N.	Error @ 100 °C	Error @ 121 °C	Error @ 130 °C
		(°C)	(°C)	(°C)
1.	1011004332	-0.07	0.12	0.08
2.	1001009966	0.19	0.32	0.28
3.	1011004203	0.23	0.17	0.09
4.	1011004509	0.36	0.16	0.19
5.	1011004211	0.51	0.53	0.48
6.	1011004106	-0.05	0.10	0.07
7.	1011004102	-0.06	0.12	0.07
8.	1011004404	0.48	0.27	0.32
9.	1011004338	0.07	0.20	0.31
10.	1011004503	0.20	0.31	0.25

Reviewed By (TAG)

Reviewed By (MEVAC)

Name

Ahmed Tarek

Signature

Date

30-06-2025



VOL-01) PERFORMANCE QUALIFICATION PROTOCOL

ANNEX 01.01

RESULTS DATA SHEETS

DATA SHEET 01.01.04

LIST OF ATTACHEMENTS



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
 Issue Date: 02/02/2025

Page 24 of 25

Datasheet Number 01.01.04.01: List of Attachments

List of Attachments		
Attachment No.	Description	Location
01.02.01.01	Maximum Solid Distribution Cycle Printouts	Annex-01.02 (Turn Over Package)
01.02.02.01	Maximum Solid Penetration Cycle Printouts	Annex-01.02 (Turn Over Package)
01.02.01.02	Raw Data (on CD)	Annex-01.02 (Turn Over Package)
01.02.01.03	Inspection Label	Annex-01.02 (Turn Over Package)
01.02.01.04	Reference Equipment Calibration Certificates	Annex-01.02 (Turn Over Package)
01.02.01.05	Inspection Certificate	Annex-01.02 (Turn Over Package)

Reviewed By (TAG)

Reviewed By (MEVAC)

Name

Ahmed Tarek

Signature



Date

30-06-2025



Jun., 2025

Rev. No.: 00

Doc. No.: PQR-V-00725

F-7.4/34
Issue Date: 02/02/2025

Page 25 of 25

FINAL REVIEW AND APPROVAL

The review of the test execution and the relevant raw data allow concluding that the performance qualification has been concluded:

Comply ☒Comply with Deviation ☐Doesn't Comply ☐

NOTES:

Approved by	<u>Salma Salah</u> Name	<u>Salma Salah</u> Signature	<u>T.O.E.</u> Title	<u>TAG</u> Company	<u>30-06-2025</u> Date
Approved by	<u>Ahmed Mohamed</u> Name	<u>Ahmed</u> Signature	<u>P.E.</u> Title	<u>TAG</u> Company	<u>30-06-2025</u> Date
Approved by	<u>Noha Essam</u> Name	<u>Noha Essam</u> Signature	<u>SR. T.O.E.</u> Title	<u>TAG</u> Company	<u>30-06-2025</u> Date
Approved by	<u>Ahmed Tarek</u> Name	<u>Ahmed Tarek</u> Signature	<u>T.M.</u> Title	<u>TAG</u> Company	<u>30-06-2025</u> Date
Approved by	 Name	 Signature	 Title	 Company	 Date
Approved by	 Name	 Signature	 Title	 Company	 Date
Approved by	 Name	 Signature	 Title	 Company	 Date

